



Clinical Assistance Guide

Aspire Supalite Boot Scooter

Aspire Supalite Boot Scooter

The Ultimate Portable Travel Scooter! Every day is a new adventure with the Boot Scooter Supalite, an exceptionally stable and manoeuvrable 4-wheel travel scooter.

Benefiting from a class-leading, ultra-lightweight design the ability to be disassembled into lightweight pieces to travel anywhere, the Boot Scooter Supalite keeps up with your active life.



FEATURES

The ultimate portable travel scooter

- An exceptionally stable and maneuverable 4-wheel travel scooter
- Class-leading, ultra-lightweight design with a max. speed up to 6.5 km/h
- Easy to transport at home, hotels, or airports
- A standard front basket provides storage

Compact seating with easy adjustments

- Compact, foldable seat with armrests
- Easy-to-grip tiller adjustment knob and handles
- Stable, front frame-mounted seat post

An ultra-lightweight pull-apart design

- Disassembles into 6 lightweight pieces for easy transport and storage
- Heaviest piece is only 10.7kg
- Up to 8kg lighter than comparable travel scooters
- Aluminium chassis construction

Rechargeable Li-Ion battery

- Removable Lithium-Ion battery back; battery can be charged on the scooter or removed for charging
- May be carried on airlines (check with the carrier)

IDEAL USER CONDITIONS

Mobility scooters are restricted to just 10km/h in most Australian States & Territories, and can only be used in pedestrian areas. As such, they're not classed as vehicles and don't need to be registered – except in Queensland, which offers free registration.

Users also don't need a licence to operate them. The only requirement in Australia is that the user has a genuine need of a mobility scooter. However, it's important that they're able to operate their mobility scooter safely. If there are concerns about physical or cognitive capacity for safe scooter use, there is benefit in an OT assessment or discussion with the GP.

Older people or people with the following abilities and support needs may benefit from using an Aspire Supalite Boot Scooter:

- Those with impaired mobility who have
 - The ability to independently transfer on/off scooter
 - Suitable vision and hearing
 - Suitable sitting balance
 - Suitable range of movement and function in their neck, shoulders, arms and hands
 - Cognitive abilities to safely manage a mobility scooter

- Pain or fatigue who would like to be able to access indoor or paved areas with level terrain
- Weigh less than the SWL of 120kgs
NB: When determining loading, it is important to take into account the user's weight plus any items that will be carried on the scooter, including, for example, accessories, O2 bottle, handbag etc
- Require a lightweight scooter that can be dismantled for travel and storage
- Would like to join family or friends at shopping centres or even on holidays. The Supalite scooter can improve access and endurance for crowded spaces and make the outing achievable

Considerations

The Aspire Supalite Boot Scooter is designed to be pulled apart and transported in the boot of a car. If the scooter user is using a van or accessible taxi (e.g. Maxi-Taxi) for transport, they are required to transfer off the scooter onto the seat of the vehicle. It is not safe to remain seated on a scooter during transport.

CLIENT OUTCOMES

The Aspire Supalite Boot Scooter is suitable for indoor use such as shopping centres or areas with level terrain. It can assist users in the following ways:

Making daily activities easier

- A decrease in mobility, whether due to age or illness, can make daily activities much harder. People may either forgo these activities entirely, or get help from others to do these things for them e.g.,
 - Assistance to do the shopping
 - Assistance to attend appointments
 - Limit or stop going out or meeting up with friends
- A mobility scooter allows the user to carry out these activities themselves, without assistance. Suddenly they are able to get around at home, travel to the shops and back, and even get to appointments, without hassle.

Increased independence and reduced carer support

- Mobility scooters allow people who may not have previously been able to travel far without assistance, to gain a significant level of independence
- Scooter use can enable independent community access from home in communities with suitable and level terrain e.g., Retirement Villages or from home to appointments/shops using a Maxi-Taxi
- The small turning circle of the Aspire Supalite Scooter enables entry to stores and businesses that have suitable access and spaces within their facility
- Use of this scooter can enable the user to attend appointments or shop independently from the carer if they choose. Carers and the scooter users can go in different directions to meet their own needs for the outing or enjoy the outing together

Community access and ability to travel

- A major benefit of a mobility scooter is that they can be taken on public transport and other accessible vehicles, such as Wheelchair Accessible Taxis. This allows the user to travel further, while being able to take their mobility scooter with them.
NB: Scooter users need to be able to transfer off the scooter and onto the taxi seat for transport.
- Mobility scooters can also be taken on holidays. Most airlines will allow users to check their scooter in as oversize luggage and the Lithium-Ion battery taken as hand luggage, allowing users greater access and endurance when on vacation
- The Supalite mobility scooter can even be left in the boot of the car so that it is available when you are out and about. Only the battery pack needs to be removed for charging

Improved health and wellbeing, choice and control

- Use of a mobility scooter can help to combat social isolation and health implications of the elderly or someone living with a disability. An activity that was previously too daunting can now seem achievable, enabling the scooter user to join family and friends in accessing the community and to lengthen the time that they can participate in an activity
- Scooter users have the ability to determine their own direction when accessing the shops or community. This can contribute to increased dignity and freedom of choice
- Reduced reliance on carers to attend appointments to assist the user with their mobility
- Improved mood resulting from social engagement, increased independence and reduced pain and fatigue

Falls prevention

- Scooter use can increase safety for those who are at risk of falling due to pain or impaired mobility

Reduced pain and fatigue

- Scooter use can support those who develop pain or fatigue with walking distances or standing e.g., walking around shopping centres; managing crowds or queueing

SPECIFICATIONS

Overall Size (L x W x H) mm	Seat Width mm	Wheel Size inch	Max Speed km/h	Ground Clearance mm	Turning Radius mm	Travel Range km	Battery Size V	Product Weight kg	Heaviest Component kg	Safe Working Load kg
1022 x 480 x 880	440	8/8	6.5	32	1385	18.5	1 x 24V 12Ah (Li-Ion)	34	10.7	120

WARRANTY

2 YEARS WARRANTY
Frame/Motors/Axles/Brakes

6 MONTHS WARRANTY
Battery

COMPLIANCE

AS3695.2:2013



MAINTENANCE

Your mobility scooter will need to be serviced annually to keep it running well. Aidacare offer mobility scooter servicing. Contact 1300 133 120 to find your local service centre.

Charging your mobility scooter:

- To get the most out of your mobility scooter, ensure the batteries are fully charged before each outing
 - Turn the scooter power switch to the OFF position
 - Connect the charger to the charging port on the scooter and plug the cord into a power outlet
 - The red/orange LED light indicates that the mobility scooter is charging
 - The light will change to green when the scooter's battery is fully charged
 - Once the green light appears, disconnect charger from the charging port and then from the power outlet
- The approximate charging time for the mobility scooter battery is 8-10 hours. We suggest charging your scooter overnight to make sure it's ready for you to take out the next day.

LINKS

- [Mobility HQ: Buyer's Guide to Mobility Scooters](#)
- [Travelling with batteries and power packs | Civil Aviation Safety Authority](#)

AIDACARE KNOWLEDGE HUB

Aidacare Website:

- [Aspire Boot Scooter Supalite Travel Scooter | Aidacare](#)

PDF Resources:

- [Aspire Supalite Boot Scooter Flyer](#)
- [Aspire Supalite Boot Scooter Clinical Guide](#)
- [Aspire Mobility Scooter Selection Guide](#)

Video Library:

- [Aspire Supalite Boot Scooter](#)