



# FLEX SEAT WALKER



- **Lightweight design** Easy to manoeuvre and transport, encouraging independence in daily activities.
- **Available in 6" and 8" wheel options** Offers flexibility to choose a frame suited to your environment and mobility needs.
- **Durable aluminium frame** Provides stability and reliability for everyday use.

- **Height-adjustable seat** Accommodates different user heights for safe sit to stand transfers.
- **Padded seat** Adds comfort for sitting breaks during daily activities.
- **Adjustable and lockable brake handles** Improve safety when walking and sitting.

**1** YEAR WARRANTY



FLEX	Overall Length mm/ in	Overall Width mm/ in	Overall Height mm/ in	Seat Depth mm/ in	Seat Width mm/ in	Seat Height mm/ in	Backrest Height mm/ in	Folded Width mm/ in	Product Weight kg/ lbs	Castor Size in
<b>6" Wheels</b> WAF709660CH Champagne WAF709660PU Lilac	730 / 28.7	570 / 22.4	730 - 830 / 28.7 - 32.6	300 / 11.8	360 / 14.1	430 - 520 / 16.9 - 20.4	710 / 27.9	290 / 11.4	6.3 / 13.8	6
<b>8" Wheels</b> WAF709670BK Black WAF709670BL Blue	820 / 32.2	585 / 23	780 - 970 / 30.7 - 38.1	300 / 11.8	360 / 14.1	475 - 570 / 18.7 - 22.4	760 / 29.2	390 / 15.3	7.3 / 16	8

ASFY9716181224