



RAPHAEL X QUATTRO ACTION

4 MOTOR LIFT RECLINE CHAIR

Designed for enhanced comfort and customisable support. The Raphael X delivers exceptional comfort with features to suit a diverse range of users. The adjustable backrest includes zipper access for personalised postural adjustment and lateral support, providing improved upright sitting. It's plush design compliments any modern home environment.



Lateral Backrest
Standard



Classic Backrest
Optional



2 YEAR
WARRANTY



SWL
130kg



MAXI SWL
180kg

Colour range



Golden Grain
Microsuede



Marina Grey
Microsuede



Nordic
Stone Vinyl



Custom
Upholstery

Functions



Back
Recline



Leg
Raise



Power
Lumbar



Power
Headrest

Supporting Comfort, Safety, and Independence in Daily Living

The Raphael X range is available in five sizes—Petite, Small, Medium, Large, and Maxi. Each size is designed to provide the right balance of comfort, postural support, and ease of use, ensuring a well-fitted seating experience for a wide range of users.



Adjustable Backrest for Personalised Support

Backrest cushions can be repositioned to accommodate individual postural needs, promoting better alignment and comfort.



Quattro Action Mechanism

Independently adjusts the backrest, leg rest, lumbar support, and headrest to provide tailored comfort for various activities, including reclining, relaxing, and seated support.



Backlit Hand Control with USB Charging

Easy-to-use, illuminated controls improve visibility in low light, while the integrated USB port allows convenient device charging without needing to leave the chair.

- **Built-in Lateral Supports:** Helps maintain an upright sitting posture, improving comfort and stability
- **Integrated Footrest Extension:** Helps with fatigue and improves lower limb comfort*
- **Adjustable Seat Depth:** Offers 25mm of depth adjustment, allowing taller users to achieve a better-seated position
- **Protective Covers for High-Wear Areas:** Removable covers for the headrest and armrests to keep your chair looking new
- **Refined, Contemporary Design:** Larger, plush armrests provide extra comfort while maintaining a streamlined look that suits home environment
- **All-day Battery Backup:** A lithium-ion battery ensures continued operation during power outages, supporting safety and independence
- **Premium Upholstery Options:** Available in modern finishes, including microsuede and leather-look vinyl, with local customisation available by request
- **Classic Backrest Option:** Available as an additional accessory for users who prefer a more traditional support profile



*Custom Order Form

Code	Colour	Size	Seat Height	Seat Depth	Seat Width	Seat to Top of Back	Overall Height	Overall Depth	Overall Width	Distance from Wall	SWL
CHP227640	Golden Grain	Petite*	450mm	450-475mm	470mm	630mm	1060mm	840mm	820mm	420mm	130kg
CHP227641	Marina Grey										
CHP227642	Nordic Stone										
CHP227650	Golden Grain	Small	450mm	480-505mm	500mm	630mm	1060mm	870mm	850mm	400mm	130kg
CHP227651	Marina Grey										
CHP227652	Nordic Stone										
CHP227660	Golden Grain	Medium	470mm	500-525mm	525mm	680mm	1100mm	900mm	875mm	400mm	130kg
CHP227661	Marina Grey										
CHP227662	Nordic Stone										
CHP227670	Golden Grain	Large	500mm	530-555mm	550mm	720mm	1100mm	910mm	900mm	420mm	130kg
CHP227671	Marina Grey										
CHP227672	Nordic Stone										
CHP227680	Golden Grain	Maxi	510mm	540-565mm	650mm	675mm	1120mm	910mm	1020mm	420mm	180kg
CHP227681	Marina Grey										
CHP227682	Nordic Stone										

*The Petite model is designed without a leg extension, catering to the needs of shorter users.