

SWING AWAY FOOTREST MOBILE SHOWER COMMUNE



Self Propelled Model



Attendant Propelled Model



View our Aspire Shower
Commode Accessory Flyer

- **Height Adjustable, Swing Away Footrests** for optimal positioning and ease of transfers
- **Angle Adjustable Backrest** from 8° to 18° for improved seating comfort
- **Ultra Low Seat Height** of 590mm which can be reduced to 490mm when customised
- **Toilet Clearance** improved over toilet clearance and 8° backrest angle ensures patient comfort
- **Integrated Access Trough** improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- **Open Armrests Style** ensures wider sitting surface and smaller base footprint for improved accessibility
- **Safety Arms** prevent falls and swing back out of the way for transfers
- **Brake Access** 4 x **TENTE** total locking castors provide brake access from all angles
- **Attendant Propelled** model available

Commode is registered in the ARTG AUST R ARTG 286695

5 YEAR WARRANTY
Corrosion

1 YEAR WARRANTY
Components



Code	Description	Height	Width	Seat Width	Backrest Height	Armrest Height	Anterior Seat Height	Posterior Seat Height
BTC065900	460mm Attendant Propelled	895 – 995mm	535mm	460mm	305mm	225mm	590mm	560mm
BTC066100	460mm Self Propelled	895 – 995mm	665mm	460mm	305mm	225mm	590mm	560mm
BTC066060	530mm Attendant Propelled	995mm	605mm	530mm	305mm	225mm	590mm	560mm
BTC066065	530mm Self Propelled	995mm	735mm	530mm	305mm	225mm	590mm	560mm
BTC061300	600mm Attendant Propelled	995mm	675mm	600mm	305mm	225mm	590mm	560mm
BTC061305	600mm Self Propelled	995mm	805mm	600mm	305mm	225mm	590mm	560mm