



# VIDA WHEELCHAIR

The Aspire VIDA is an extremely lightweight wheelchair that can be lifted in and out of a vehicle easily by carers. The VIDA range comes in various fun, bold colours and is ideal for quick trips to the shops, appointments, and catching up with friends. The ergonomic backrest and seat shape mimic body contours for comfort and postural support. The high safe working load (135kg / 21st) caters to a wide range of users.



Minimum Chair Weight from **7.5kg**



**UPHOLSTERY**  
Mustard/Black  
**FRAME**  
Grey  
**SEAT WIDTH**  
400mm



**UPHOLSTERY**  
Teal/Black  
**FRAME**  
Grey  
**SEAT WIDTH**  
450mm



**UPHOLSTERY**  
Magenta/Black  
**FRAME**  
Grey  
**SEAT WIDTH**  
450mm



**UPHOLSTERY**  
Olive/Black  
**FRAME**  
Grey  
**SEAT WIDTH**  
500mm



**5** YEAR WARRANTY  
Frame

**1** YEAR WARRANTY  
Components

Available as attendant propelled or self-propelled  
 Three seat width sizes: 400mm, 450mm, and 500mm, each with its own unique colour  
 The 400mm seat width is available in mustard  
 The 450mm seat width is available in teal and magenta  
 The 500mm seat width is available in olive  
 Lightweight frame design ensures easy lifting in and out of cars, ideal for carers  
 Contoured and padded seating improves comfort and postural support  
 Attendant operated hand brakes with parking lock  
 Simple and intuitive folding for easy transportation  
 Removable leg rests and folding back for compact storage  
 Rear storage pocket for essentials  
 Quick release wheels  
 Anti-tippers included for added stability  
 Lap belt available (WCA848300)



### Attendant Propelled (AP) & Self Propelled (SP) Vida Specifications

400mm Mustard

450mm Teal or Magenta

500mm Olive



MWS449874 - AP  
 MWS449876 - SP

MWS449878 - AP  
 MWS449880 - SP

MWS449882 - AP  
 MWS449884 - SP

MWS449886 - AP  
 MWS449888 - SP



View our  
 Aspire Wheelchair  
 Accessories Range

**VIDA**  
 WHEELCHAIR

| Overall Length | Overall Width | Overall Height | Seat Depth | Seat Width | Anterior Seat Height | Posterior Seat Height | Backrest Height | Armrest Height | Min Chair Weight | Rear Wheel |
|----------------|---------------|----------------|------------|------------|----------------------|-----------------------|-----------------|----------------|------------------|------------|
| mm/in          | mm/in         | mm/in          | mm/in      | mm/in      | mm/in                | mm/in                 | mm/in           | mm/in          | kg/lbs           | in         |

#### Attendant Propelled

|           |        |          |          |          |        |          |        |        |       |          |    |
|-----------|--------|----------|----------|----------|--------|----------|--------|--------|-------|----------|----|
| MWS449874 | 995/39 | 600/23.5 | 930/36.5 | 420/16.5 | 400/16 | 500/19.5 | 460/18 | 430/17 | 200/8 | 7.5/16.5 | 13 |
| MWS449878 | 995/39 | 650/25.5 | 930/36.5 | 420/16.5 | 450/18 | 500/19.5 | 460/18 | 430/17 | 200/8 | 7.5/16.5 | 13 |
| MWS449882 | 995/39 | 650/25.5 | 930/36.5 | 420/16.5 | 450/18 | 500/19.5 | 460/18 | 430/17 | 200/8 | 7.5/16.5 | 13 |
| MWS449886 | 995/39 | 680/27   | 930/36.5 | 420/16.5 | 500/20 | 500/19.5 | 460/18 | 430/17 | 200/8 | 8/17.6   | 13 |

#### Self Propelled

|           |         |        |          |          |        |          |        |        |       |          |    |
|-----------|---------|--------|----------|----------|--------|----------|--------|--------|-------|----------|----|
| MWS449876 | 1040/41 | 630/25 | 950/37.4 | 420/16.5 | 400/16 | 500/19.5 | 460/18 | 430/17 | 200/8 | 7.5/16.5 | 22 |
| MWS449880 | 1040/41 | 690/27 | 950/37.4 | 420/16.5 | 450/18 | 500/19.5 | 460/18 | 430/17 | 200/8 | 7.5/16.5 | 22 |
| MWS449884 | 1040/41 | 690/27 | 950/37.4 | 420/16.5 | 450/18 | 500/19.5 | 460/18 | 430/17 | 200/8 | 7.5/16.5 | 22 |
| MWS449888 | 1040/41 | 710/28 | 950/37.4 | 420/16.5 | 500/20 | 500/19.5 | 460/18 | 430/17 | 200/8 | 8/17.6   | 22 |

Note: Range of motion dimensions are accurate to a +/-10mm tolerance

Australia  
 1300 133 120 | aspirecare.com.au

United Kingdom  
 01622 541235 | aidacare.co.uk