

## SHOWER STOOL WIDE



- Clip on seat attachment for easy assembly and cleaning
- Height adjustable legs with easy push button height adjustment for changing needs
- · All frames individually treated for corrosion resistance
- Seat with anterior tilt that helps facilitate standing from sitting and provides greater comfort for stiff or painful hips
- Raised armrests for assistance in transferring on and off the stool
- Durable non-slip feet provide excellent grip and stability in wet conditions and a stainless steel washer enhances durability
- Additional height increments on one leg to accommodate uneven surfaces
- Comfortable wide seat with drainage holes for increased support and hygiene
- Reinforced wide frame to provide stability for larger users
- Tested to AS/NZS3973.7:2009







SHOWER STOOL wide	Overall Height	Overall Width	Overall Depth	Seat Height	Seat Width	Seat Depth	Between Arms	Product Weight
	mm	mm	mm	mm	mm	mm	mm	kg
STEEL BTS119010	575 - 765	570	500	425 - 615	420	365	525	3.8
ALUMINIUM BTS119015	575 - 765	570	500	425 - 615	420	365	525	3.1



OVA CONTRACTED



AS / NZS Approved

